


Making the Invisible Visible

Dynamic Nuclear Polarization (DNP) experiments transfer the high polarization of electron spins to nuclear spins, driven by microwave irradiation of unpaired electron spins. Bruker DNP-NMR spectrometers are designed specifically for extended solid-state NMR experiments, delivering unsurpassed sensitivity and stability for exciting new applications in biological solids, material science and pharmaceuticals.



Gyrotron and Klystron Microwave Sources

Bruker's custom-designed gyrotron microwave sources and low-temperature MAS (LTMAS) probes have a proven record of performance with 38 installed systems to date. We also offer a 263 GHz Extended Interaction Klystron (EIK) with 5 W output power and high frequency/power stability. The 263 GHz klystron provides 80-100% of the gyrotron DNP efficiency (depending on the sample) with a lower purchase price, footprint and facility requirements.

Magnetic Field	¹ H NMR Frequency	EPR/μwave Frequency	Microwave Source
9.4 T	400 MHz	263 GHz	Klystron
14.1 T	600 MHz	395 GHz	Gyrotron
18.8 T	800 MHz	527 GHz	Gyrotron
21.1 T	900 MHz	593 GHz	Gyrotron





Fig. 1: Polymeric carbon nitride materials doped with 15 mM AMUPol in aqueous solvent yield excellent DNP signal enhancements of 71 with ~5 W of microwave power, allowing ¹⁵N-¹³C correlation spectra to be acquired rapidly even at natural isotopic abundance.

1.3 and 1.9 mm MAS DNP Probes

The DNP LTMAS probes operate in the 100-200 K temperature range with cold insert/eject capabilities. They are offered with HCN, HX, HXY (with variety of X/Y combinations) or low gamma RF configuration and the following rotor sizes to cover a range of applications at 400 to 900 MHz:

- 3.2 mm: 15 kHz MAS @ 100 K
- 1.9 mm: 24 kHz MAS @ 100 K
- 1.3 mm: 40 kHz MAS @ 100 K

The 1.3 mm probe has optimized microwave coupling into the sample, providing high DNP efficiency with fast MAS.




Fig. 2: EM simulations for 1.3 mm optimized DNP stator and waveguide at 263 GHz.

Fast MAS + DNP

As in conventional solid-state NMR, fast magic angle rotation enables superior decoupling of dipolar interactions, often provides favorable relaxation properties and can even allow for ¹H-detection. Combined with large signal enhancements from DNP, this has opened up new applications in structural biology and materials science.

1.3 mm DNP probes provide ideal sensitivity for limited-quantity samples. As in the case of a ¹³C,¹⁵N-specifically amino acid labeled dihydrofolate reductase, 4 μL of sample at a concentration of 0.65 mM were sufficient to acquire a ¹H-¹³C heteronuclear correlation spectrum in just 26 minutes.




Fig. 3: ¹H-¹³C HETCOR spectrum of specifically-¹³C,¹⁵N-labeled dihydrofolate reductase at a concentration of only 0.65 mM (2.6 nanomoles of protein in rotor), enhanced with 20 mM TOTAPOL in a 3:7 v/v glycerol-d₈/D₂O buffer, 40 kHz MAS.

Summary

- Turn-key solution for DNP-enhanced solids NMR experiments at high field.
- Unique high power CW gyrotron microwave sources at 263, 395, 527, 593 GHz.
- Klystron microwave source option at 263 GHz for increased DNP accessibility.
- Low-temperature (100 K) MAS probe technology with built-in waveguide and cold spinning gas supply.
- High DNP signal enhancements on wide range of samples.